

# Point Endo

## Product Comparison



### Point Endo

The Point Endo is a two joint ratcheting prosthetic finger specifically designed to be integrated into a cosmetic silicone covering. It is designed for **light to medium duty applications** and has a narrower profile and smooth exterior with no grip surfaces or contours.



### Point Digit

The Point Digit is a three joint ratcheting prosthetic finger to be used without any cosmetic covering. It is designed for **heavy-duty applications** and features a wider profile with a contoured exterior and integrated grip surfaces.



---

Both the Point Endo and the Point Digit share a few key features that aid in durability and function. These include their ability to be unilaterally operated, anatomically flexed around the MCP joint, and to withstand loads up to 150 lbs.



### Choosing the Right Product

Given their shared features, both the Point Endo and the Point Digit offer similar functionality. There are key differences though, which are listed below.

#### Point Endo

- Light/medium-duty applications
- Compatibility with custom cosmetic cover
- Lower price point
- Only available outside of the U.S.\*
- Fewer locking positions
- Smaller range of motion

#### Point Digit

- Heavy-duty applications
- Not compatible with a custom cosmetic cover
- Higher price point
- Available in the U.S. and internationally
- More locking positions
- Larger range of motion

\* Please reference the Point Endo Distributors document to search for product availability