



# Point Partial Order Form

PO#:

Needed On/Before

Fill out a sizing form for **EACH** hand being fitted. Return the form via e-mail to [sales@pointdesignsllc.com](mailto:sales@pointdesignsllc.com).

**This form must be filled out completely before the order can be filled.**

## Clinic

Clinic Name

Clinician's Name

Email

Phone

## Shipping

Street Address / P.O. Box

City, State

ZIP code

## Patient

Patient Name

Dominant Hand?

Left

Right

Affected Hand?

Left

Right

Affected Finger(s)?

2 (index)

3 (middle)

4 (ring)

5 (little)

Please list the functional expectations for your device

List the 5 most frequent manual tasks to be performed. (i.e., typing, playing piano, cutting food, etc.)

Patient Height

Patient Weight

Patient Age

Today's Date

Date of Amputation



point designs

### Sizing

*please complete each step below for proper sizing:*

1. Measure the distances from the PIP joint centers to the fingertips on the intact hand (where applicable) and record in the table to the right.
2. Consider socket build out and any residual limb distal to the PIP joint when choosing sizes. In general, round down to the nearest size.
3. Choose a size from the table below for each desired prosthetic partial digit.
4. Choose a finger surface finish option and a mounting kit from the table below (45 mm, 50 mm, 55 mm).

**Overall Length (mm)**  
*PIP Joint Center to Distal Fingertip*

Index	<input type="text"/>	mm
Middle	<input type="text"/>	mm
Ring	<input type="text"/>	mm
Little	<input type="text"/>	mm

*Please include additional material if you feel it would be useful (e.g., photographs, additional drawings and measurements, photographs with rulers, scans, etc.)*

### Digit Size Table

<u>Length</u>	<u>Part Number</u>	<u>Index</u>	<u>Middle</u>	<u>Ring</u>	<u>Little</u>
45 mm	PNTPL2-045-G	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
50 mm	PNTPL2-050-G	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
55 mm	PNTPL2-055-G	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Finger surface finish option

<input type="checkbox"/> Gunmetal Titanium
--

### Mounting Kit

<u>Size</u>	<u>Part Number</u>	<u>Quantity</u>
45 mm	PNTPL2MK-045	
50 mm	PNTPL2MK-050	
55 mm	PNTPL2MK-055	

**Are you willing to conduct a quickDASH outcome measure prior to the fitting and after the fitting?**

A \$5 Amazon gift card will be provided after submission of the pre-fitting quickDASH, and a \$15 Amazon gift card will be provided after submission of the post-fitting quickDASH.

Yes  No

*The above information is true to the best of my knowledge. By signing the form, I am accepting responsibility for the sizing. If the sizing is incorrect, resulting in an ill-fitting device, I understand I may incur additional charges.*

Clinician Signature \_\_\_\_\_ Date \_\_\_\_\_