



The


POINT PIVOT

by point designs



point designs

POINT PIVOT



Rotate Clockwise

Push Digit
Down



Rotate Counterclockwise



Compact Configuration



Extended Configuration



point designs

The Point Pivot is a rotatable mounting system that provides an additional degree of freedom to the Point Thumb. As a replacement to the standard mounting system for the Point Thumb, the Point Pivot provides allows rotation of the flexion/extension plane of the Point Thumb, allowing for a greater number of achievable hand grasps.

The Point Pivot comes in two configurations: Compact and Extended. Compact is designed for integration with new sockets. Extended is designed for integration with existing sockets where a Point Thumb was previously mounted.