

The

POINT PIVOT+


by point designs



point designs

www.pointdesignsllc.com | sales@pointdesignsllc.com


POINT PIVOT+



point designs

The Point Pivot+ is a dynamic mounting system that provides additional degrees of freedom to the Point Thumb. As a replacement to the standard mounting system for the Point Thumb, the Point Pivot+ provides many functional benefits to the patient. It allows the Point Thumb to be moved from an oppositional position to a lateral position. Additionally, it enables rotation of the flexion/extension plane of the Point Thumb, allowing for a greater number of achievable hand grasps.

Adduction/Abduction



Rotation

